# HOW TO CREATE THE IDEAL GELMA FORMULATION FOR SPECIFIC APPLICATIONS A story of temperature and photo-crosslinking

#### Jan-Philip Zegwaart<sup>1,2</sup>, Thomas Van Gansbeke<sup>1</sup>, Catarina Ferreira da Silva<sup>1</sup>, Riccardo Levato<sup>2</sup> & Jos Olijve<sup>1</sup>

1. Rousselot Biomedical, Meulestedekaai 81, 9000, Gent, Belgium. 2. Department of Orthopaedics, University Medical Center Utrecht, Utrecht University, 3584 CX Utrecht, The Netherlands

## INTRODUCTION


Gelatins are widely used in the biomedical and tissue engineering fields, due to their bio-compatibility, biodegradability and cell-interactivity. Gelatin can be chemically crosslinked to originate hydrogel constructs that are stable at body temperature. However, to achieve **predictable hydrogel strength** several important factors need to be considered. Starting with **consistent molecular weight (MW)** and **degree of modification (DoM)** levels. Next, the resin formulation needs to be on point; the interplay of (bio)polymer concentration, salt concentration, photo-initiator (PI) type and concentration, are curically important. Because gelatin is thermo-responsive, **temperature control in the pre-phase and during hydrogel production**, is essential. The interplay is diverse, but mappable and insightful decisions can be made early in the research phase.

### PURPOSE

The herein presented study has as aim to provide some elucidation and insight into the **effects of temperature**, in hydrogel production, in the pre-phase and during the **photo-crosslinking phase** of methacrylated gelatins (**GeIMA**).

## MATERIALS AND METHODS

Four GeIMA types 90p40; 90p60; 160p60; 160p80 (90 kDa at DoM of 40 and 60 %, and 160 kDa at DoM of 60 and 80 %), were dissolved at 5 % w/v in 1x PBS. The LAP photoinitiator was used at a 0.1% w/v concentration in the GeIMA resin formulations. The crosslinking kinetics and the hydrogel strengths of the GeIMAs were studied at various temperatures (45°C; 37°C; 21°C & 4°C) using photo-rheology (A 25 mm plate-plate setup with a quartz bottom plate; 1 mm gap; MCR 302e, Anton Paar, Belgium). Oscillation measurement (frequency: 1 Hz; Shear strain: 1%) with curing: the temperature was set to 45°C; 37°C; 21°C or 4°C (tolerance 0.1°C for 10 min to allow temperature equilibration). A parafin oil corona was applied to prevent sample drying. The G' and G" were measured for 2 minutes. Then, G' and G" were measured during the curing of the sample for 45 minutes at 11 mW cm<sup>-2</sup> at a wavelength of 365 nm. After, G' and G" were measured for an additional 7 minutes.



 4°C
 affected by the temperature during the curing process.<sup>10</sup>

 4°C
 11 Larry L. Hench and Ion K. West, The sol-gen process, Chamila Reviews 1990 90 (1), 33-72, DOI: 10.1021/cr00099a003;

 12] S. Panja, B. Dietrich, D. J. Adams, Controlling Symersis of Hydrogel Stype Organic Salts, Angew. Chem. Int. Ed. 2022, 61, e202115021;

 4°C
 regeneration, Bioactive Materials 15 (2022) 82-102, https://doi.org/10.1016/j.bioactmat.2021.12.005

With consistent GelMA molar mass and degree of modification, hydrogel strength becomes predictable. However, the applied protocols greatly affect hydrogel outcomes, as demonstrated herein through the temperatures applied in the preparation and the photo-curing phases. In short, the interplay between the various components of a GelMA resin is diverse, but it is mappable and hence insightful decisions can be made, greatly improving reproducibility and results for biomedical scientists.

#### Acknowledgement

0

37°C

This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 964497

21°C





